**STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL
REINFORCEMENT IMPLANTS**

en	Instructions for use	Page	2
fr	Notices d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
pt	Nota de instruções	Pagina	12
ελ	Οδηγίες χρήσης	σελίδα	14
da	Brugsvejledning	Side	16
fi	Käyttöohjeet	Sivu	18
hu	Felhasználási útmutató	Pagina	20
nl	Gebruiksaanwijzing	Pagina	22
ro	Instrucțiuni de utilizare	Pagină	24
tu	Kullanım Kilavuzu	Sayfa	26
pl	Instrukcja użycia	Strona	28
cs	Pokyny k použití	Strana	30
ru	Инструкция по применению	страница	32
sk	Návod k použití	Strana	34
sl	Navodila za uporabo	Stran	36
sv	Bruksanvisning	Sida	38
no	Bruksanvisning	Side	40
uk	Інструкція із застосування	сторінка	42
zn	使用说明	页	44
ar	تعليمات الاستخدام	صفحة	46

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This release is the last update of the instructions of use and replace the previous edition

PROMESH® SURG ST
STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT
IMPLANTS
SINGLE USE STERILE PRODUCT

DESCRIPTION

PROMESH® SURG ST meshes are non-absorbable parietal reinforcement implants.

IMPLANTED MATERIALS

Mesh: polypropylene (knitted monofilament), blue thread (depending on references): monofilament polypropylene non absorbable, colorant [phthalocyaninato(2-)] copper.

Origin neither human nor animal – Non-absorbable.

INDICATIONS

Repair and parietal reinforcement for inguinal and femoral hernias and eventrations.

PERFORMANCE

Meshes are particularly adapted either to laparoscopy surgery, or to conventional surgery.

The PROMESH® SURG ST is designed for extraperitoneal implantation only.

The mesh is biocompatible and non-absorbable, causing a fibrosis reaction which takes over the reinforcement after six months of implantation. It has the advantage of having a shape memory, easy to use dimensions, a very high resistance to suturing, and great flexibility.

The blue thread on reshaped meshes is the anatomic landmark for the pubic bone

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant could generate possible undesirable side effects:

- | | | | |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma | - Foreign Body Reaction |
| - Infection | - Fistula formation | - Mesh migration | - Irritation near by organ |
| - Recurrence | - Erosion | - Mesh deformation | - obstruction |
| - Inflammation | - Seroma | - Allergic reaction | |

USAGE PRECAUTIONS

PROMESH® SURG ST meshes are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN SURGERY does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

When the surgeon considers it necessary to do so, it is recommended to cut the implant outside the surgical site to avoid any risk of releasing particles into the patient.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the **PROMESH® SURG ST** mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

Store in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN SURGERY representative or distributor.

EXPLANTATION AND ELIMINATION OF DEVICES

Explanations and handling should be carried out in accordance with ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN SURGERY on request. It is important to note that any implant that has not been cleaned and disinfected before dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUEST AND COMPLAINTS

In accordance with its quality policy, COUSIN SURGERY is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN SURGERY as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN SURGERY without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN SURGERY and its distributors.

Caution : Federal law (USA) restricts this device to sale,distribution and use by or on the order of a physician.

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