# LIGAMENT VERTEBRAL LIGANOVE

---

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	5
de	<a href="#">Gebrauchsanweisung</a>	Seite	8
it	<a href="#">Istruzioni per l'uso</a>	Pagina	11
es	<a href="#">Instrucciones de uso</a>	Pagina	14
pt	<a href="#">Nota de instruções</a>	Página	17
kr	<a href="#">사용설명서</a>	페이지	20

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Made in France

**NOT087 / 260421**  
Version du 21/04/2026

This release is the last update of the instructions of use and replace the previous edition

**DEVICE FOR THE INTER-SPINOUS SPACE WITH SPINOUS SUPPORT  
LIGANOVE VERTEBRAL LIGAMENT**

Reinforcement ligaments for spinal surgery Articular ligament  
STERILE PRODUCT FOR SINGLE USE ONLY

**Description**

The LIGANOVE vertebral ligament is a medical device for use in spinal surgery.

**Materials**

Polyethylene terephthalate (polyester) – Barium sulphate – Polyvinyl chloride (see packaging) - No human or animal origin – Non-absorbable.

**Indications**

The LIGANOVE vertebral ligament is recommended for use in thoracolumbar spine (from T1 to L5 levels) in the following cases:

- Spinal surgery –
- Artificial Vertebral ligament –
- Dorsolumbar instability

**Performances**

The main function of the LIGANOVE vertebral ligament is to assist the posterior articular ligament during flexion.

**Contraindications**

Do not implant the device in the following cases:

- Allergy to one of the components –
- Infection at the surgical site –
- Children that have not finished growing
- Pregnant women
- Surgery at more than 2 levels with LIGANOVE vertebral ligament

**Undesirable secondary effects**

As with all implantable medical devices, this implant may cause undesirable secondary effects; these may include:

- Infection
- Device breakage (Rupture of the ligament)
- Uncrimping of the needle
- Device loosening, Permanent elongation of the ligament
- Device migration
- Seroma/swelling around the device
- Inflammatory reactions, allergy, foreign body reaction
- Persistent or worsening pain
- Spinous process fracture
- Recurrence with revision surgery or conversion to fusion
- Wrong level implantation, inappropriate or inaccurate positioning
- Removal or revision surgery, Surgical revision with intersomatic arthrodesis in the operated area
- Hypersensitivity reaction

Complications unlikely to be related to the device but related to surgery or to the patient's conditions:

- Infection
- Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis
- Loss of neurological function
- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain
- Urinary retention or loss of bladder control or other types of urological system compromise
- Fracture, Bone loss or decrease in bone density

- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union. Mal-union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise
- Development of respiratory problems

### **Precautions for use**

The LIGANOVE vertebral ligament is supplied sterile (sterilized with ethylene oxide or gamma radiation). Prior to use, check the integrity of the device and its packaging (blister packs or peel-open pouches). If the device or its packaging is damaged in any way, do not use the device. Do not use if the device is past its expiry date.

COUSIN BIOTECH cannot provide any guarantees or recommendations regarding the use of a specific brand of fixation device. The properties of these devices are subject to modifications by the manufacturers, over which COUSIN BIOTECH has no control. Avoid the use of staples or screws with sharp edges.

This device should only be implanted by a qualified surgeon trained in the use of the product.

LIGANOVE is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

### **IMPORTANT: DO NOT REUSE - DO NOT RE-STERILIZE**

As indicated on the product labelling, this implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness).

### **Precautions for storing the prosthesis: Store in a dry place, away from sunlight, at room temperature, and in its original packaging.**

Should you require any additional information regarding the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **MRI SAFETY:**

The implants are composed of non-ferromagnetic materials and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori, they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of the implant, has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded a conditional compatibility for MRI scan up to 1.5 Tesla.

As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, the patient who has this implant must be advised to warn the health professionals concerned (radiologists and radiology operators) in as far as possible about the presence of this implant before these exams.

### **"EXPLANTATION AND ELIMINATION OF DEVICES"**

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

## **“INFORMATION REQUEST AND COMPLAINTS”**

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country, and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

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