# PROMESH® SURG UMB

## STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

| | | | | |
|---|---|---|---|---|
| en | Instructions for use | Page | 2 | *CE marking date: 09/05/2025* |
| fr | Notices d'instructions | Page | 6 | *Date de marquage CE: 09/05/2025* |
| de | Gebrauchsanweisung | Seite | 10 | *Datum der CE-Kennzeichnung: 09/05/2025* |
| it | Istruzioni per l'uso | Pagina | 14 | *Data di marcatura CE: 09/05/2025* |
| es | Instrucciones de uso | Página | 18 | *Fecha de marcado CE: 09/05/2025* |

Manufactured by **COUSIN BIOTECH**
Allée des Roses
59117 Wervicq-Sud – France
Tél : +33 (0) 3 20 14 41 20
Fax : +33 (0) 3 20 14 40 13
http://www.cousin-biotech.com

Distributed by **Peters Surgical**
Immeuble AURELIUM
1 cours de l'île Seguin
92100 BOULOGNE BILLANCOURT-FRANCE
Tel : +33 (0) 1 48 10 62 62
Fax : +33 (0) 1 48 91 22 99
https://www.peters-surgical.com

Made in France

CE 0476

**NOT312_261026**
Version du 26/10/2026

℞ only Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition
**Product under regulation 2017/745 (EU).**

**STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS**

**PROMESH® SURG UMB**

EXPANDABLE PARIETAL PROSTHESIS WITH INFLATABLE BALLOON FOR VISCERAL SURGERY

**SINGLE USE STERILE PRODUCT**

## MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

## IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The **PROMESH® SURG UMB** prosthesis consists of:
- A parietal prosthesis made of an ePTFE and polypropylene plug, two or four polyethylene terephthalate with biocompatible colorant (D&C Green #6) and ePTFE fixation threads crimped with stainless steel needles (see on package).
- A balloon which allows the deployment of the parietal prosthesis (except for PSUMBIN05D).

## INDICATIONS

Repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non absorbable support material.

## CONTRAINDICATIONS

Do not use in the following cases:
- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

## UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

| | | | |
|---|---|---|---|
| - Discomfort/Pain | - Recurrence | - Adhesion formation | -Chronic Pain |
| - Erosions | - Irritation nearby organ | - Infection | -Organ injury |
| - Inflammation | - Fistula formation | - Seroma | |
| - Mesh deformation | - Hematoma | - Mesh migration | |
| - Allergic reaction | - Foreign body reaction | - Ileus /Obstruction | |

## TARGET POPULATION

Adult with ventral hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.

## IMPLANTED MATERIALS

Polypropylene, Polyethylene terephthalate, ePTFE, Biocompatible colorant (D&C Green #6)

| Reference | Weight of ePTFE (g) | Weight of Polypropylene (g) | Weight of PET (g) |
|---|---|---|---|
| PSUMBIN05D | 0,01 g | 0,57 g | 0,19 g |
| PSUMBIN07D | 0,02 g | 1,13 g | 0,19 g |
| PSUMBIN09D | 0,04 g | 2,03 g | 0,19 g |

## MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY :

Stainless steel, Poly vinyl chloride.
**Origin neither human nor animal – Non-absorbable**.

## PERFORMANCES

The **PROMESH® SURG UMB** prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for ventral, umbilical and incisional hernias.
Products, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow a optimum and rapid integration and colonisation.
If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation.

## LINK TO THE SAFETY AND CLINICAL PERFORMANCE SUMMARY :

https://www.ifu-oem.com/d-62c425f69bacd_summary-of-safety-and-clinical-performance-of-sterile-intraperitoneal-non-resorbable-parietal-reinforcement-implants.pdf

## MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene, polyethylene terephthalate and non metallic component. Consequently, **PROMESH® SURG UMB** is considered MR safe.

## PRECAUTIONS FOR USE

**PROMESH® SURG UMB** prosthesis are delivered sterile (ethylene oxide sterilization).
Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.
Do not use if the device is out of date.
COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.
This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
**IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE**

### THIS BALLOON MUST IMPERATIVELY BE WITHDRAWN

As specified on the product's labeling, the **PROMESH® SURG UMB** mesh is for single use only. It cannot be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

## INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

### STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

### SURGICAL TECHNIQUE

| SURGICAL PROCEDURE (except for PSUMBIN05D which comes without balloon) | | | | | |
|---|---|---|---|---|---|
| **1-Installation:**<br>- The prosthesis equipped with its balloon crosses completely the hernial defect,<br>-its positioning becomes then entirely retro parietal (intraperitoneal) | | **2-Deployment:**<br>- The expansion and the spreading of the prosthesis are obtained by inflating the balloon,<br>- Parietal fixation is carried out by two sutures thanks to the pre-positioned threads | | **3- Prosthesis in place:**<br>-The prosthesis is fixed,<br>-The balloon must imperatively be withdrawn | |

## Description of the prosthesis:

This prosthesis is in a shape of a flat plug with 2 fixation threads crimped with needles (4 needles) for the product references PSUMBIN05D and PSUMBIN07D, and with 4 fixation threads crimped with needles (8 needles) for the PSUMBIN09D product reference.

The aim of its innovating expansion balloon concept is to make dissection easier and to secure the unfolding of the prosthesis.
This balloon is located inside the plug; it is filled-up with 30cm$^3$ of air. It acts as a pneumatic Jack ensuring complete crease-free unfolding of the prosthesis.
Its double-sided design combines a smooth ePTFE side which reduces viscera adhesions with a knitted monofilament polypropylene side which allows better colonisation of the peritoneum.

**Sizes of the prosthesis:**
The medical device is a sterile kit.
The product range is composed of:
PSUMBIN 5cm diameter: for small hernias and ventral hernias on trocart site
PSUMBIN 7cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias < 2cm
PSUMBIN 9cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias that measure between 2 and 5cm

**Installation of the patient :**
Patient lying on his back
Both legs kept in the body's axis
Each arm placed on an arm wrester in a 90° angle from the body's axis

**1-Incision :**
The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia.

**2-Dissection and resection of the peritoneal sac :**

The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place.
Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand.

**3-Positioning and fixating the prosthesis :**
**PROMESH® SURG UMB** is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way.
The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice.
The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall.
Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

The balloon is deflated and removed.

**4-Closure of the incision :**
Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

## INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.
Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

## IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

## EXPLANTATION AND DISPOSAL OF DEVICES

This device is intended for permanent parietal reinforcement and to remain permanently implanted. Routine removal or replacement is not required. In the event of certain complications, a secondary surgical procedure for partial or complete mesh explantation may be considered, if clinically indicated. Complete removal of all implant material may not always be technically feasible. The decision to explant the device lies with the responsible surgeon and shall be based on the individual clinical situation and current medical practice.
Mesh explantation is a technically demanding reoperation and should only be performed by a surgeon who is appropriately qualified, trained, and experienced in mesh removal techniques. As with any reoperation, mesh explantation carries significant risks, including but not limited to bleeding, infection, injury to adjacent anatomical structures (e.g. nerves, vessels, bowel, bladder, or other organs), anaesthesia-related complications, recurrence or persistence of the original condition, worsening of symptoms, chronic pain, and the need for additional surgical interventions. Patients should be informed that mesh explantation may resolve, improve, partially improve, fail to resolve or worsen the reported symptoms or adverse effects.
If a hernia mesh implant is explanted and retrieval data are to be collected, the device should be retrieved and handled in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants – Part 1: Retrieval and handling". Explanted devices that are not returned for analysis must be disposed of in accordance with applicable local regulations for infectious and hazardous clinical waste.


**IMPLANT CARD INFORMATION LEAFLET**

## How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

**Front side of the implant card**

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

**Back side of the implant card**

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

| Information to identify the device | Traceability information |
|---|---|
| • Indication that the product is a medical device<br>• Description of the device family<br>• Trade name of the device<br>• Reference of the device<br>• Device is MR safe<br>• Denomination name of the device in sales countries | • UDI number<br>• Batch number<br>• Name and address of the legal manufacturer |

## Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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