PROMESH® SURG DOME

STERILE SEMI-RESORBABLE REINFORCEMENT PARIETAL IMPLANTS

en	Instructions for use	Page	2	*MDR CE marking date: 30.09.2025*
fr	Notices d'instructions	Page	6	*Date de marquage CE MDR : 30.09.2025*

Manufactured by **COUSIN BIOTECH**
ALLEE DES ROSES
59117 Wervicq-Sud – France
Tél : +33 (0) 3 20 14 41 20
Fax : +33 (0) 3 20 14 40 13
http://www.cousin-biotech.com

Made in France

Distributed by **Peters Surgical**
Immeuble AURELIUM
1 Cours de l'Ile Seguin
92100 Boulogne-Billancourt - FRANCE
Tel : +33 (0) 1 48 10 62 62
Fax : +33 (0) 1 48 91 22 99
https://www.peters-surgical.com

CE 0297

NOT303_260929
Version du 29/09/2026

℞ only Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

Instructions for use *en*

PROMESH® SURG DOME

SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT
SINGLE USE STERILE PRODUCT

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The PROMESH® SURG DOME device consists of a semi-resorbable dome and a semi-resorbable reinforcement mesh for anterior and posterior parietal reinforcement of the groin.

INTENDED USE

Intended to be used as parietal reinforcement implant.

INDICATIONS

Repair and surgical reinforcement of inguinal and femoral hernias.

CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Recurrence
- Hematoma
- Foreign body reaction
- Infection
- Adhesion
- Mesh migration
- Irritation nearby organ
- Inflammation
- Erosion/Extrusion
- Mesh deformation
- Ileus / Obstruction
- Fistula
- Seroma/lymphocele
- Allergic reaction
- Organ injury
- Chronic Pain
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes

TARGET POPULATION

Adult male or female suffering from groin hernias with an indication of surgical repair with extraperitoneal implantation. These include:

- Inguinal hernias
- Femoral hernias

IMPLANTED MATERIALS

Polypropylene (PP) (non resorbable) - Poly-L-Lactic Acid (PLLA) (resorbable)

		Weight PP (g)	Weight PLLA (g)
PSABDO2PCL	**Dome**	0,03	0,26
	Flat mesh	0,21	0,64
PSABDO3PCL	**Dome**	0,05	0,41
	Flat mesh	0,21	0,64
PSABDO4PCL	**Dome**	0,08	0,60
	Flat mesh	0,21	0,64
PSABDO2PCS	**Dome**	0,03	0,26
	Flat mesh	0,12	0,37
PSABDO3PCS	**Dome**	0,05	0,41
	Flat mesh	0,12	0,37
PSABDO4PCS	**Dome**	0,08	0,60
	Flat mesh	0,12	0,37

Origin neither human nor animal origin. Semi-resorbable

PERFORMANCES

The device PROMESH® SURG DOME contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, PROMESH® SURG DOMEis considered MR safe.

PRECAUTIONS FOR USE

The dome should always be installed with an anterior reinforcement.
The dome should not be cut.
PROMESH® SURG DOME devices must always be separated from abdominal cavity by peritoneum
PROMESH® SURG DOME devices are delivered sterile (ethylene oxide sterilisation).
When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.
Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.
Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.
Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.
Prevent the mesh from any contact with objects that may alter their surface.
This device must be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the PROMESH® SURG DOME implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section “Implant card information leaflet” for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss**. The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

PROMESH® SURG DOME Surgical Technique (1)

Indirect Hernia

❶ In indirect inguinal hernia, the sac is dissected up to the internal inguinal ring

❷ The PROMESH® SURG DOME is placed in the defect maintaining the sac deeply. The PROMESH® SURG DOME is then fixed by 1 or 4 stitches with non-resorbable suture to the resistant edge of the internal ring

❸ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in the Lichtenstein technique.

PROMESH® SURG DOME Surgical Technique (2)

Direct Hernia

❶ In direct hernia the sac is dissected up to the transversalis fascia.

❷ The PROMESH® SURG DOME is placed in the opening whilst holding the sac deeply. The PROMESH® SURG DOME is then fixed by 7 or 8 stitches with non-resorbable suture to the resistant edge of the transversalis fascia.

❸ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in Lichtenstein's technique.

EXPLANTATION AND DISPOSAL OF DEVICES

This device is intended for permanent parietal reinforcement and to remain permanently implanted. Routine removal or replacement is not required. In the event of certain complications, a secondary surgical procedure for partial or complete mesh explantation may be considered, if clinically indicated. Complete removal of all implant material may not always be technically feasible. The decision to explant the device lies with the responsible surgeon and shall be based on the individual clinical situation and current medical practice.

Mesh explantation is a technically demanding reoperation and should only be performed by a surgeon who is appropriately qualified, trained, and experienced in mesh removal techniques. As with any reoperation, mesh explantation carries significant risks, including but not limited to bleeding, infection, injury to adjacent anatomical structures (e.g. nerves, vessels, bowel, bladder, or other organs), anaesthesia related complications, recurrence or persistence of the original condition, worsening of symptoms, chronic pain, and the need for additional surgical interventions. Patients should be informed that mesh explantation may resolve, improve, partially improve, fail to resolve or worsen the reported symptoms or adverse effects.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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